

**Description of this Thermometer**

- ① ON/OFF button
- ② Display
- ③ Battery compartment cover
- ④ Measuring sensor / measuring tip
- ⑤ Cleaning and disinfecting area (thermometer probe only)

EN

Antimicrobial Copper Properties

Surfaces made or covered by special copper alloys, have strong antimicrobial properties against a wide variety of microorganisms*. Copper alloys emit antimicrobial copper ions Cu+ that whilst in contact with microbes and bacteria rupture their cellular membranes, thus destroying these microorganisms. This activity reduces the microbial flora on the coated area and – due to the "halo phenomenon" – simultaneously causes a drastic reduction in pathogens on the remaining body of the thermometer. Thermometers containing parts with antimicrobial copper alloys drastically reduce microbial flora, minimizing the dispersion of contagious microorganisms and thus providing high safety to the end user.

*Efstathia A. Panos "The Role of Antimicrobial Copper Surfaces in Reducing Healthcare associated Infections", European Infectious Disease, Volume 5, Issue 2, Autumn 2011.

Important Safety Instructions

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- This device is only to be used for measuring human body temperature through oral, rectal or axillary. Do not attempt to take temperatures at other sites, such as in the ear, as it may result in false readings and may lead to injury.
- Do not use the device if you think it is damaged or notice anything unusual.
- We recommend cleaning this device according to the cleaning instructions before first use for personal hygiene.
- The minimum measurement time until the beep is heard must be maintained without exception!
- Consider that different measurement locations may require continued measuring even after the beep, see section "Measuring methods / Normal body temperature".

For safety reasons (risk of perforation) rectal measurement in children younger than 3 years must be performed only by trained health care personnel (professional user). Use another measurement method instead. For rectal fever measurement in small children younger than 3 years, thermometers with a flexible tip are available.

• Do not attempt rectal measurements on persons with rectal disorders. Doing so may aggravate or worsen the disorder.

• Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

• Do not use this device close to strong electromagnetic fields such as mobile phones or mobile installations. Keep a minimum distance of 3 cm from the device when using this device.

• Protect the device from impact and drops.

• Avoid ambient temperatures above 60 °C. NEVER boil this device!

• Use only the commercial disinfectants listed in the section "Cleaning and Disinfecting".

• We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact your local Microlife Service to arrange the test.

WARNING: The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Read the instructions carefully before using this device.

Type BF applied part

Turning on the Thermometer

To turn on the thermometer, press the ON/OFF button ①, a short beep signals "thermometer On". A display test is performed. All segments should be displayed.

Then at an ambient temperature of less than 32 °C, an "L" and a flashing "C" appear at the display field ②. The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by "ERR" on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

Using the Thermometer

Before use, keep this thermometer out of any physical contact, for at least 2 hours. This utilizes antimicrobial copper properties, as described in chapter "Antimicrobial Copper Properties". The use of antimicrobial copper is a supplement to and not a substitute for standard infectious control practices. Users must continue to follow all current infection control and cleaning practices. We recommend cleaning the thermometer as described in chapter "Cleaning and Disinfecting".

The preferred measurement location is the axilla. When taking a measurement, the measurement location is continuous on the display and a "C" symbol flashes. If the beep is heard 10 times in a row, "C" is no longer flashing, this means that the measured increase in temperature is less than 0.1 °C in 16 seconds and the thermometer can now be read.

To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button ⑤. Otherwise the thermometer will automatically turn off after about 10 minutes.

Storage of Measured Values

If the ON/OFF button ① is pressed for more than 3 seconds when turning on the thermometer, the automatically stored maximum temperature during the last measurement will be displayed. At the same time, a "M" for memory will appear on the display. About 2 seconds after the button is released, the temperature value disappears and the thermometer is ready for measurement.

Measuring methods / Normal body temperature

In the armpit (axillary) | 34.7 - 37.3 °C

Attention: For the prevention of rectal perforation in children (younger than 3 years), we recommend using another measuring method, or the use of a thermometer with flexible tip.

Carefully insert the measuring sensor ④ of the thermometer 2 to 3 cm into the anal aperture.

The use of a probe cover and the use of a lubricant is recommended. If you are unsure of this measurement method, you should consult a professional for guidance/training.

Approx. measurement time: 1 minute!

In the anus (rectal) | 36.6 - 38.0 °C

Attention: For the prevention of rectal perforation in children (younger than 3 years), we recommend using another measuring method, or the use of a thermometer with flexible tip.

Carefully insert the measuring sensor ④ of the thermometer 2 to 3 cm into the anal aperture.

The use of a probe cover and the use of a lubricant is recommended. If you are unsure of this measurement method, you should consult a professional for guidance/training.

Approx. measurement time: 1 minute!

Cleaning and Disinfecting

For disinfection in home use environment, use a 70% isopropyl alcohol swab, or a cotton tissue moistened with 70% isopropyl alcohol to wipe surface pollutants off the thermometer probe (note: consider the application and safety instruction of the disinfectant manufacturer).

Always start wiping from the end of the thermometer probe (approx. at the middle of the thermometer) towards the thermometer tip. Afterwards the entire thermometer probe (see number ⑤ in the drawing)

should be immersed in 70% isopropyl alcohol for at least 5 minutes (max. 24 hours). After immersion, let the disinfectant dry off for 1 minute before next use. Avoid immersing or wiping the display to protect it from fading.

The thermometer is not intended for professional use.

Battery Replacement

When the "▼" symbol (upside-down triangle) appears in the display, the battery is flat and needs replacing. To replace the battery remove the battery compartment cover ③ from the thermometer. Insert the new battery with the + at the top. Make sure you have a battery of the same type to hand. Batteries can be purchased at any electrical store.

Technical Specifications

Type: Maximum thermometer
Measurement range: 32.0 °C to 43.9 °C
Temp. < 32.0 °C: display "L" for low (too low)
Temp. > 43.9 °C: display "H" for high (too high)

Measurement accuracy: ± 0.1 °C, 34 °C - 42 °C
± 0.2 °C, 32.0 - 33.9 °C and 42.1 - 43.9 °C

Operating conditions: 10 - 40 °C, 15-95 % relative maximum humidity

Storage conditions: -25 - +60 °C, 15-95 % relative maximum humidity

Battery: LR41 (1.5V) / SR41 (1.55V)
Battery lifetime: approx. 4500 measurements (using a new battery)

IP Class: IP22

Reference to standards: EN 12470-3, medical thermometers; ASTM E1112, IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11

Expected service life: 5 years or 10000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC.

Technical alterations reserved.

Guarantee

This device is covered by a 5 year guarantee from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee.

• Consider that different measurement locations may require continued measuring even after the beep, see section "Measuring methods / Normal body temperature".

• For safety reasons (risk of perforation) rectal measurement in children younger than 3 years must be performed only by trained health care personnel (professional user). Use another measurement method instead. For rectal fever measurement in small children younger than 3 years, thermometers with a flexible tip are available.

• Do not attempt rectal measurements on persons with rectal disorders. Doing so may aggravate or worsen the disorder.

• Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

• Do not use this device close to strong electromagnetic fields such as mobile phones or mobile installations. Keep a minimum distance of 3 cm from the device when using this device.

• Protect the device from impact and drops.

• Avoid ambient temperatures above 60 °C. NEVER boil this device!

• Use only the commercial disinfectants listed in the section "Cleaning and Disinfecting".

• We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact your local Microlife Service to arrange the test.

WARNING: The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Read the instructions carefully before using this device.

Type BF applied part

Guarantee

This device is covered by a 5 year guarantee from the date of purchase.

During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

The following items are excluded from the guarantee:

• Transport costs and risks of transport.

• Damage caused by incorrect application or non-compliance with the instructions for use.

• Damage caused by accident or misuse.

• Packaging/storage material and instructions for use.

• Regular checks and maintenance (calibration).

• Accessories and wearing parts. Battery.

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service through our website: www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

Opis termometru

• Przycisk ON/OFF (wl/wyl.)

• Wyświetlacz

• Zatyczka pojemnika na baterię

• Czujnik pomiarowy / końcówka pomiarowa

• Obszar czyszczenia i dezynfekcji (tylko sonda termometru)

• Czujnik pomiarowy / końcówka pomiarowa

• Wykrywanie przekroczenia temperatury

• Wykrywanie przekroczenia czasu

</div



Popis teploméru

- ① Tlačítko ON/OFF (ZAP/VYP)
- ② Displej
- ③ Kryt prostoru pro baterie
- ④ Měřicí senzor / měřicí hrot
- ⑤ Čisticí a dezinfiční prostor (pouze sonda teploméru)

Antimikrobiální vlastnosti médi

Povrchy vyrobené ze speciálního slitin médi nebo nimi pokryté, mají silné antimikrobiální vlastnosti proti širokemu spektru mikroorganismů. Slitiny médi uvoľňujú antimikrobiální mědi Cu, které v kontaktu s mikromi a bakteřemi způsobují prasknutí buňkových membrán, čím týmikroorganizmy zmrzlí. V důsledku toho dochází k redukcii patogenů na ostatních částech teploméru. Teplomer, které obsahují dily s antimikrobiálními slitinami radikálně redukují mikrobiální flóru, minimalizují riziko napakly mikroorganismů, a tím zajišťují vysokou bezpečnost koncového užívateľa.

E. Saito A. Panos / The Role of Antimicrobial Copper Surfaces in Reducing Healthcare associated Infections, European Infectious Disease, Volume 5, Issue 2, Autumn 2011.

Důležité bezpečnostní pokyny

• Postupujte podle návodu k použití. Tento návod obsahuje důležité informace o chodbu a bezpečnosti tohto prístroja. Před používáním prístroje si důkladně prečtěte tento návod a uschovejte jej po další použití.

• Tento prístroj je možné použíti pouze pro měření lidské telesné teploty orálně, rectálně alebo axilárně (v podpaží). Nepoužíte se měřit teplotu na jiných miestach, kde by došlo k uchlu, protože to mohlo vést k chybám výsledku měření a může to způsobit zranění.

• Nepoužívejte tento prístroj, pokud si myslíte, že je poškozený nebo má výraznější nádor.

• Pro osobní hygienu doporučujeme před prvním použitím přistrop odstranit a nechat ho na čistě.

• Vždy je nutno dodržet minimální dobnu měření, dokud se neozve písmotu.

Vezměte na vědomí, že různá místa měření mohou vyzádovat počítání měření i po zaznamenání zvukového signálu, viz kapitola "Metody měření / Normální telesná teplota".

• Bezpečnostních úvodů (riziko perforace reakta) musí rektální měření v délce 3 měsíců 3 roky provádět pouze skouzené zdravotnický personál (profesionální uživatel). Vzdušný prostředek. Vzdušný prostředek.

• Práva a ochrana osobního soukromí.

• Používejte tento prístroj, aby ho nechal souborníkem.

• Používejte tento prístroj, aby ho nechal souborníkem.